


PRODUCT DATA SHEET

Revision 1	35mm Confocal Cell Culture Dish	
Validfor Item-No.:	1050001/1051001/1050000/1051000/ 1050011/1051011/1050010/1051010	

1.	Description / Specification	
1.1	Description	35 mm Confocal Cell Culture dish, glass bottom, sterile 1050001 35 mm Confocal Cell Culture Dish , hole diameter 15 mm, TC treated, sterile 1051001 35 mm Confocal Cell Culture Dish, hole diameter 20 mm, TC treated, sterile 1050000 35 mm Confocal Cell Culture Dish, hole diameter 15 mm, TC treated, sterile 1051000 35 mm Confocal Cell Culture Dish, hole diameter 20 mm, TC treated, sterile 1050011 35 mm Confocal Cell Culture Dish, black ,hole diameter 15 mm, TC treated, sterile 1051011 35 mm Confocal Cell Culture Dish, black , hole diameter 20 mm, TC treated, sterile 1050010 35 mm Confocal Cell Culture Dish, black hole diameter 15 mm, TC treated, sterile 1051010 35 mm Confocal Cell Culture Dish, black hole diameter 20 mm, TC treated, sterile
1.2	Dimensions	<u>Dish</u> : see Customer Drawing
1.3	Volume	Max. volume: 4 ml Working volume: 3 ml Growth area: 15mm (1.766 cm ²) and 20mm (3.14 cm ²)
1.4	Material / Resin	<u>Dish</u> : PS (Polystyrene) <u>Lid</u> : PS (Polystyrene) <u>Bottom</u> : Borosilicate Glasses The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	<u>Dish</u> : clear or black (1050011/1051011/1050010/1051010) <u>Lid</u> : clear <u>Bottom</u> : clear
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. Printed on Packaging

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-
3. Packaging		
3.1	Pieces / Bag	10/20
3.2	Pieces / Box	200
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download
4. Other Information		
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.