


PRODUCT DATA SHEET

Revision 1	ELISA Microplate, 384 Well, Flat-Clear, Undetachable	
Validfor Item-No.:	2072010/2072000/2072011/2072001/2073011/2073001	

1.	Description / Specification	
1.1	Description	ELISA Microplate, 384 Well, PS, Flat-Clear, Undetachable 2072010 ELISA plate, 384 well, flat-bottom clear, Undetachable, medium binding force 2072000 ELISA plate, 384 well, flat-bottom clear, Undetachable, high binding force 2072011 ELISA plate, 384 well, flat-bottom clear, Undetachable, high binding force, without lid 2072001 ELISA plate, 384 well, flat-bottom clear, Undetachable, middle binding force, without lid 2073011 ELISA plate, 384 well, flat-bottom clear, Undetachable, middle binding force, with lid 2073001 ELISA plate, 384 well, flat-bottom clear, Undetachable, high binding force, with lid
1.2	Dimensions	Plate: see Customer Drawing
1.3	Volume	Max. volume: 131 µl Working volume: 15 -110 µl
1.4	Material / Resin	<u>Plate</u> : PS (Polystyrene) <u>Lid</u> : PS (Polystyrene) The materials for manufacturing are USP Class VI
1.5	Colour	<u>Plate</u> : clear <u>Lid</u> : clear <u>Bottom</u> : clear
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for immunology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the Packaging

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	NO
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	2072010/2072000	2072011/2072001/2073011/2073001
3.1	Pieces / Bag	10	1
3.2	Pieces / Box	200	100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 07.05.2024	Date 14.05.2024	Date 31.05.2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.