


PRODUCT DATA SHEET

Revision 1	ELISA Microplate, 96 Well, Flat-Black, Undetachable	
Validfor Item-No.:	2070210/2070200/2070211/2070201	

1.	Description / Specification	
1.1	Description	ELISA Microplate, 96 Well, PS, Flat-Black, Undetachable 2070210 ELISA plate, 96 well, flat-black, Undetachable, medium binding force 2070200 ELISA plate, 96 well, flat-black, Undetachable, high binding force 2070211 ELISA plate, 96 well, flat-black, Undetachable, middle binding force, without lid 2070201 ELISA plate, 96 well, flat-black, Undetachable, high binding force, without lid
1.2	Dimensions	Plate: see Customer Drawing
1.3	Volume	Max. volume: 388 µl Working volume: 20-350 µl
1.4	Material / Resin	Plate: PS (Polystyrene) Lid: PS (Polystyrene) The materials for manufacturing are USP Class VI
1.5	Colour	Plate : black Lid: clear Bottom: black
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for immunology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the Packaging

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	NO
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	2070211/2070201	2070210/2070200
3.1	Pieces / Bag	1	10
3.2	Pieces / Box	100	200
3.3	Lot-No.	DDD YY XXX (day, year, product No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 07.05.2024	Date 14.05.2024	Date 31.05.2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.