


PRODUCT DATA SHEET

| | | |
|--------------------|------------------------------|---|
| Revision 1 | 96 Well ELISA Microplate Lid |  |
| Validfor Item-No.: | 2073000 | |

| 1. | Description / Specification | |
|-----|-----------------------------|---|
| 1.1 | Description | 96 Well ELISA Microplate Lid, sterile 2073000 ELISA plate cover, 96 well, clear, flat-bottom, sterile |
| 1.2 | Dimensions | Lid: see Customer Drawing |
| 1.3 | Volume | Max. Volume: - Working volume: - |
| 1.4 | Material / Resin | Lid: PS (Polystyrene) The materials for manufacturing are USP Class VI |
| 1.5 | Colour | Lid: clear |
| 1.6 | Sterilization | No |
| 1.7 | Quality Control | - Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Intended Use | General laboratory product for immunology to be used by qualified personnel in a laboratory environment. |
| 1.9 | Other Information | - For single use only - Expiry date and Lot-No. printed on the Packaging |

| 2. | Features | |
|-----|--------------------------|--|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic |
| 2.2 | Temperature range | For application: -20°C to +60°C |
| 2.3 | Autoclavability | NO |
| 2.4 | Centrifugation, max. RCF | N/A |
| 2.5 | Shelf life | 3 years |
| 2.6 | Other Information | - |

| 3. | Packaging | |
|-----|-------------------|-------------------------------------|
| 3.1 | Pieces / Bag | 25 |
| 3.2 | Pieces / Box | 750 |
| 3.3 | Lot-No. | DDD YY XXX (day, year, product No.) |
| 3.4 | Other Information | Certificate of Quality to download |

| 4. | Other Information | |
|-----|---------------------|--|
| 4.1 | Laboratory use only | |

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved. |
|--------------------|---------------------|--------------------|---------------------|---|
| Revision 1 | Date 07.05.2024 | Date 14.05.2024 | Date 31.05.2024 | |
| Date 01.06.2024 | Name Chen Anjing | Name Zhang Heng | Name Wang Kedong | |