


# PRODUCT DATA SHEET

Revision 1	UV Analyser Plate	
Validfor Item-No.:	2110001/2111001	

1.	Description / Specification	
1.1	Description	UV Analyser Plate 2110001 Ultraviolet analysis plate, 96 well, independent blister 2111001 Ultraviolet analysis plate, 384 well, independent blister
1.2	Dimensions	Plate : see Customer Drawing
1.3	Volume	Max volume: 131 µl /388 µl Working volume: 15 -110 µl /20-350 µl
1.4	Material / Resin	Plate : PS (Polystyrene) The materials for manufacturing are USP Class VI
1.5	Colour	<u>Plate:</u> clear
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for molecular biology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the Packaging

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 07.05.2024	Date 14.05.2024	Date 31.05.2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.